

## **INTRODUCING PORTATANK®**

Portatank® is the ideal solution for non-invasively identifying levels of water and clean liquids in large industrial tanks

An innovative handheld tool to accurately identify the level of contents in large containers.

- **Type** Portable Ultrasonic Liquid Level Indicator
- Function Designed for inspecting the liquid levels of tanks with widths between 0.5
   15 metres
- **Part Number** 3141184-TANK

Quick Testing
Very clear LED results
for liquid indication

Versatile
Applicable for a wide range of tanks & liquids

Non-Invasive
Testing outside tanks
through the outer wall



## **ADVANTAGES VS. OTHER LEVEL METHODS**



QUICK: easy to use, instant liquid level detection



ACCURATE: to +/-10mm with very clear LIQUID / NO LIQUID results



SAFE: Non-invasive, non-intrusive — suitable for testing hazardous liquids



EFFICIENT: Handheld, portable device - no installation time or costs required



NO RISK: Measures levels of fluids, particularly hazardous liquids, where the tank cannot be opened for safety reasons and a float gauge cannot be installed



NO DOWNTIME: liquid level measurement can be taken when the tank systems are "live"



ADAPTABLE: can work with liquid levels that fluctuate over time



FLEXIBLE: suitable for a wide range of cylinders & liquids

### **APPLICATIONS**

Quickly, safely and easily identify the liquid level in containers between 0.5-15 metres diameter range.

**Tank Diameter Range:** 0.5 metres – 15 metres

Tank Wall Material: Most common types such as steel, plastic, aluminium etc.

Liquid types: Most standard industrial liquids. E.g. water and oils, if you are unsure

please ask our experts

- For tanks that already have an in-built level system, our portable non-invasive technology allows operators to verify liquid levels in tanks with existing in-built level indicators, as some of these can be from an old technology that can be inaccurate such as float level or sight-glass methods.
- For tanks without any form of liquid level indication, the Portatank® can detect the liquid levels non-invasively through the outer walls of the tank.
- Portatank® requires no installation for either application

Fire Safety Engineering Processing Pharmaceuticals

Oil & Gas

Water Treatment Facilities


Chemicals Petrochemicals

Get in touch with your requirements to find out how we can help your

application!

## **HOW TO TEST: 3 SIMPLE STEPS**

After unit has been set up for the application.



- 1. Apply ultrasonic gel to the tank wall, connect the sensor to the unit and turn the unit on
- 2. Place the sensor on the tank wall starting from the upper section, then move the sensor down in steps.
- 3. The red LED indicator will light up when liquid level is encountered

# **HOW DOES PORTATANK® WORK**

Clear results: The red LED indicator will light up when liquid level is encountered. Clear LIQUID / NO LIQUID indication.

Working principle: Ultrasonic pulseecho measurement. An ultrasonic signal is transmitted through the walls, across the tank's contents and its reflection is measured to determine the presence or absence of liquid.

**Troubleshooting made easy:** by connecting the included oscilloscope



### ADVANCED TROUBLESHOOTING



### **Handheld Oscilloscope Specification**

• Screen size: 2.4 inch

• Bandwidth: 30 MHz

• Sampling rate: 200 MS/s

Vertical sensitivity: 50mV – 100V

• Time range: 50S – 15nS

• Trigger mode: Single / Normal / Auto

Maximum test voltage: 1X = 60V, 10X = 800V

• Voltage precision: +/- 2%

• Battery: 1200 mAh, rechargeable via micro

• USB

• Standby time: 3 hours

#### WHY USE OSCILLOSCOPE:

- 1. Allows you to SEE the ultrasound signal travelling across the tank
- 2. Helps you SELECT the optimum parameters for liquid level detection specific to your tank, enhances accuracy and reliability
- 3. Assists you in UNDERSTANDING the underlying issues associated with any unsuccessful measurement

### WHEN IS THE OSCILLOSCOPE REQUIRED:

- 1. Setting up the Portatank® with your application for the first time: To select the optimum parameters for your tank for enhanced accuracy and reliability
- 2. Investigate the underlying issues if an unsuccessful measurement is made:

  To view the ultrasound signal and identify the appropriate solution


DISCLAIMER: THE OSCILLOSCOPE IS NOT REQUIRED DURING NORMAL OPERATION OF THE PORTATANK- ONCE THE UNIT IS SET UP TO WORK FOR YOUR APPLICATION.

# **TECHNICAL SPECIFICATIONS**

Dimensions	Main Unit: (H)168mm x(W)79mm x(D)44mm
Weight	Main Unit: 520 grams (with batteries)
Accuracy	+/- 10mm or +/-5mm (sensor dependent)
Verifiable Agents	Water and most clean liquids (get in touch for more details)
Sensor options	200 kHz, 500 kHz, 1 MHz, 2 MHz (depending on application) connected by BNC connector
Power	4 x 1.5V AA batteries
Display	Liquid Level Detection: Visual LED indication
Operating Temperature	Ambient temperature range: -20°C to +70°C
Warranty	Main Unit: 3 Years Warranty Sensor: 1 Year Warranty
Package Contents	<ul> <li>1 Main Unit</li> <li>1 Sensor</li> <li>1 Handheld Oscilloscope (+ BNC cable &amp; Micro USB Cable)</li> <li>1 Robust Carrying Case</li> <li>1 Calibration Certificate</li> <li>2 Operating Instructions</li> <li>1 Ultrasonic Gel</li> </ul>
Certificate	CE, Manufactured in UK with ISO 9001:2015 registered company

### **PACKAGE**

Comes in a ready to go package so that your servicing teams will be able to start testing as soon as they have read the Operating Instructions.



Holdoff Time (to control the period Gain (to increase incoming signal when the device starts to listen for strength to trigger LED indicator, the echo, please refer to manual)

please refer to manual)

### Also included but not pictured:

- Carrying Case
- Instruction Manual
- Calibration Certificate
- Ultrasonic Gel
- BNC cable (for connecting oscilloscope)
- Micro USB cable (for charging oscilloscope)


Technical Support provided free of charge for the unit's lifetime.

Every unit provided by Coltraco comes with a 3 year warranty supporting the manufacturing quality of the main unit and 1 year on sensor.

OPTIONAL PORTACARE® after sales package for Total Care maintenance.

We have local partners to support you worldwide through our global network of Partners, Distributors and ODA (Organisational Delegated Authorities).

ODAS are Service Centres which provide regional support for caliibration in UK, France, Greece, UAE, India, Singapore and USA.



Coltraco®, Coltraco North America®, Portamarine®, Portalevel®,Permalevel®, Portagauge®, Portasonic®, Portamonitor™, Portasteele™, Portascanner®, Permascanner®, SafesiteTM, SafesiteTM are trademarks or registered trademarks of Coltraco Limited, UK. DuPont ™, FM-200®, FE-25™, FE-13™, and FE-241™ are trademarks or registered trademarks of E.I. du Pont de Nemours and Company and its affiliates. Novec™ 1230 is a trademark owned by 3M.

**COLTRACO ULTRASONICS** 

**REVIEWS.**io

3 REVIEWS.io